**USER INSTRUCTIONS**

# INDICAID™

## COVID-19/FLU A&B Rapid Antigen Test

**Suitable for ages 2+ years**  
**Must be ages 14+ to use kit unsupervised**  
**For in vitro diagnostic use only**  
 Read all instructions carefully before performing the test. Failure to follow the instructions may result in inaccurate test results.

**INTENDED USE**

INDICAID™ COVID-19/FLU A&B RAPID ANTIGEN TEST is an in vitro diagnostic test for determining the presence of SARS-CoV-2, influenza A (flu A) and/or influenza B (flu B) antigens in anterior nasal swab samples. This test is intended for self-collected anterior nasal swab from individuals over 14 years old or adult collected anterior nasal swab from individuals aged 2 years or older with symptoms of respiratory infections within the first 6 days of symptom onset.

Results are for the identification of SARS-CoV-2, flu A and/or flu B antigens. Antigens are generally detectable in anterior nasal swabs during the acute phase of infection and are parts of the virus that can serve as markers for disease exposure.

**CONTAINED IN THIS BOX**

P0098	P0099
1x Individually-Wrapped Test Device, Buffer Solution Vial, Individually-Wrapped Swab, User Instructions	12x Individually-Wrapped Test Devices, Buffer Solution Vials, Individually-Wrapped Swabs; 1x User Instructions

**Performing Your Test**

- Check the expiration date** on the outside of the product box.
- Wash your hands thoroughly for at least 20 seconds before and after testing.
- Remove the test device and swab from their packaging.
- Tilt your head back. Gently insert the swab about **2 cm** into one of your nostrils. Rub the swab against the wall of one of your nostrils at least **4 times** in a large circular path. Repeat with your other nostril **using the same swab**.
- The buffer solution vial cap is composed of two parts. Remove the entire cap. Immediately place the nasal swab into the buffer solution vial. Tilt the vial and stir the swab into the buffer solution by **twisting** the swab back and forth **20 times**. Ensure the swab tip (soft end) is fully submerged in the solution. **Press and roll** the swab tip against the inner wall of the vial to remove excess solution.
- Close the entire vial cap tightly. Remove the top half of the vial cap to expose the dropper tip.
- Place the test device on a flat surface**. Hold the vial vertically. Slowly squeeze 4 drops of the solution into the sample well (S). Inaccurate results may occur if less than **4 drops** are applied.
- Leave for **15 minutes** and read the results. Do not read after 20 minutes. Refer to the "Interpreting Your Results" section below.
- Dispose of all used test kit supplies and swab samples in a trash receptacle.

**!** Make sure the test device is placed on a horizontal (flat) surface while running the test/until the timer ends.

**LIMITATIONS OF THE TEST**

- Positive results do not rule out bacterial infection or co-infection with other viruses and the agent detected may not be the definite cause of disease.
- Negative results do not rule out SARS-CoV-2, influenza A and/or influenza B infections, especially if you have been in contact with the viruses. A follow-up PCR test should be considered to rule out infection.
- False negative result may occur if the level of antigen in the sample is below the test limit or the test procedure is not carried out properly.
- Results from this test should not be used as the sole basis to diagnose or exclude COVID-19, influenza A and/or influenza B infections. User should not take any decision of medical relevance without first consulting the healthcare professional.
- This test kit can test for current infection only and cannot tell if you have had COVID-19, influenza A and/or influenza B infection in the past.
- A negative test result may occur if virus has mutated at the specific antigen region recognized by this test.


**WARNINGS PRECAUTIONS AND SAFETY INFORMATION**

- ✓ Children aged 2 to 13 years should be tested by an adult.
- ✓ Use only the contents provided in the test kit.
- ✓ Keep test kit and kit components away from children and pets before and after use. Avoid contact with your skin, eyes, nose, or mouth.
- ✓ Leave test device sealed in its pouch until just before use. Once opened, the test device should be used within 2 hours.
- ✓ Store at 2-30°C. Do not freeze. Avoid direct sunlight.
- ✓ Wear a face mask or other face covering while collecting sample from a child or an adult.
- ✓ The reagent solution contains hazardous chemicals. If contact to the body occurs, flush with copious amount of water. If irritation persist, seek medical advice.
- ✓ Keep foreign substances and household cleaning products away from the test device during the testing process. Contact with foreign substances and household cleaning products (e.g., 1% bleach) may result in an incorrect test result.
- ✗ Do not touch the swab tip.
- ✗ Do not reuse. Test components are for single-use.
- ✗ Do not use test kit beyond its printed expiration date.
- ✗ Do not use if any of the test kit content or packaging is damaged or opened.
- ✗ Do not ingest any kit components.
- ✗ Do not use the test on children under 2 years of age.

**Get help here:** please contact Customer Service at +852 3700 8888 or email to [cs@indicaid.com](mailto:cs@indicaid.com)

**Interpreting Your Results**

- Look for lines next to the 'C' (Control), 'A', 'B', and 'CoV'. Use the information on the right to interpret the lines.
- **Look closely!** Lines can vary in color and may appear very faint!
- **Make sure there is a line next to the 'C'!** Your result is INVALID if there is no 'C' line.

INVALID TEST RESULT No 'C' Line	NEGATIVE TEST RESULT 'C' Line Only	POSITIVE TEST RESULT 'C' Line AND any 'A', 'B', or 'CoV' Line
<p><b>STOP</b> If no 'C' line is seen, the test result is <b>INVALID</b> even if you see 'A', 'B', or 'CoV' line(s).</p> <p><b>Take these next steps:</b></p> <ul style="list-style-type: none"> <li>• Collect a new nasal swab sample and repeat the test with a new INDICAID™ COVID-19/FLU A&amp;B Rapid Antigen Test.</li> <li>• If you develop COVID-19, flu A and/or flu B symptoms or your symptoms become severe, seek medical attention immediately.</li> </ul>	<p><b>Take these next steps:</b></p> <ul style="list-style-type: none"> <li>• If you develop COVID-19, flu A, or flu B symptoms or your symptoms become severe, seek medical attention immediately.</li> <li>• If you have no symptoms, you are recommended to take a PCR test to confirm your negative result if you have recent exposure to known or suspected COVID-19 case.</li> </ul>	<p>Flu A Positive    Flu B Positive    COVID-19 Positive</p> <p>Flu A &amp; COVID-19 Positive    Flu B &amp; COVID-19 Positive</p> <p><b>Look very closely!</b> Even a very faint line next to the 'A', 'B', or 'CoV' should be considered positive if there is also a 'C' line.</p> <p><b>Take these next steps:</b></p> <ul style="list-style-type: none"> <li>• Consult your healthcare provider to discuss your positive test result.</li> <li>• Self-isolate at home per Department of Health's recommendations to stop spreading the virus to others.</li> </ul>

**Understanding Your Results**

**What an invalid result means:** The test could not tell whether or not you have COVID-19, flu A, or flu B. The test needs to be repeated with a new kit and freshly collected anterior nasal swab sample.

**What a negative result means:** The proteins from SARS-CoV-2, flu A, and/or flu B were not detected in the sample. A negative result does not rule out SARS-CoV-2, flu A and/or flu B infection. There is a higher chance of false negative results with antigen tests than with laboratory-based molecular tests. This means that there is a higher chance this test will give you a negative result when you have SARS-CoV-2, flu A, and/or flu B. If you tested negative and continue to experience COVID-19, flu A, or flu B-like symptoms of fever, cough, and/or shortness of breath, you should seek follow-up care with your healthcare provider.

**What a positive result means:** The SARS-CoV-2, flu A and/or flu B virus(es) were detected in your sample. It is very likely that you have the respective infection(s) and are contagious. Please contact your doctor and follow the Department of Health's instructions regarding the positive cases arrangement. There is a very small chance that this test can give an incorrect positive result (false positive). Your healthcare provider will work with you to determine how best to care for you based on your test results along with your medical history and symptoms.

Co-infection with SARS-CoV-2, flu A and/or flu B is rare. If results are positive for more than one test lines, the sample should be retested with new test kit. Report all results obtained from this test kit to your healthcare provider.

**PHASE Scientific International Limited**  
 1/F, Building 22E, Phase 3, Hong Kong Science Park,  
 Shatin, New Territories, Hong Kong

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© 2022 PHASE Scientific International Ltd., 版權所有